

SURGICAL KITS



D I T R O N D E N T A L U S A

DRILL STOPPER TUBE



- Segment Lengths: 8, 10, 11.5 & 13 mm
- Segment Diameters: 2, 2.5, 2.8, 3.2 & 3.8 mm

FULL DRILL STOPPER KIT



FULL DRILL STOPPER KIT: LAYOUT

Marking
Drill



6.0mm

8.0mm

10mm

11.5mm


13mm

16mm

	Ø2.0mm	Ø2.5mm	Ø2.8mm	Ø3.2mm	Ø3.8mm	Ø4.5mm
L6.0mm						
L8.0mm						
L10.0mm						
L11.5mm						
L13.0mm						


DO NOT EXCEED 130°C (274°F) DURING STERILIZATION

FULL SURGICAL KIT




- Marking Drill: Diameter 1.9 mm
- Pilot Drill: Diameter 2 mm
- Twist Drills: Diameters 2.5, 2.8, 3.2, 3.8, 4.5 & 5 mm
- Countersink Drills: Standard & Wide
- Drill Extender
- Instruments: Parallel Pin, Short & Long Hand Key, Short & Long Ratchet Key, Short & Long Motor Adapter, Short & Long Hand & Ratchet Key, Short & Long Motor Insertion Tool, Depth Probe & Ratchet Wrench

FULL SURGICAL KIT: LAYOUT



DO NOT EXCEED 134° C (273°F) DURING STERILIZATION

FULL SURGICAL KIT WITH TORQUE RACHET



- Marking Drill: Diameter 1.9 mm
- Pilot Drill: Diameter 2 mm
- Twist Drills: Diameters 2.5, 2.8, 3.2, 3.8, 4.5 & 5 mm
- Countersink Drills: Standard & Wide
- Drill Extender
- Instruments: Parallel Pin, Short & Long Hand Key, Short & Long Rachet Key, Short & Long Motor Adapter, Short & Long Hand & Rachet Key, Short & Long Motor Insertion Tool, Depth Probe & Rachet Wrench
- Torque Rachet: 10-45 Ncm

MINI SURGICAL KIT



- Marking Drill: Diameter 1.9 mm
- Pilot Drill: Diameter 2 mm
- Twist Drills: Diameters 2.5, 2.8, 3.2 & 3.8 mm
- Countersink Drill: Standard
- Drill Extender
- Instruments: Short & Long Ratchet Key, Short & Long Motor Adaptor, Short & Long Hand & Ratchet Key, Short & Long Motor Insertion Tool & Ratchet Wrench

MINI SURGICAL KIT: LAYOUT



SURGICAL TUBE



- Cover
- Mid-Tray: 7 Silicon Inlets for Drills & Instruments
- Base: Dappen-Dish
- Marking Drill: Diameter 1.9 mm
- Pilot Drill: Diameter 2 mm
- Twist Drills: Diameters 2.8, 3.2 & 3.8 mm
- Instruments: Long Ratchet Key, Long Hand & Ratchet Key & Short Ratchet Wench

ORDER NOW



D I T R O N  D E N T A L U S A